Serial break

Reset the system. The system serial communication settings return to the default: 9600 bauds, 8 data bits, no parity, 1 stop bit, no hardware handshaking.

REPLY

RESET<CRC16><CR>

<CR> is '\r'

COMM

Set serial communication settings for system

MODE

All modes

SYNTAX

COMM:<BaudRate><DataBits><Parity><StopBits><HardwareHandshaking><CRC16>
<CR>

BaudRate: 1 byte. 0 for 9600 bps, and 5 for 115200 bps.

DataBits: 1 byte. 0 for 8 bits by default, and 1 for 7 bits.

Parity: 1 byte. 0 for none. StopBits: 1 byte. 0 for 1 bit.

HardwareHandshaking: 1 byte. 0 for off.

REPLY

OKAY<CRC16><CR>

ERROR<ErrorCode><CRC16><CR>

ErrorCode: 2 bytes.

OTHER

Wait at least 100ms

<u>INIT</u>

Initialize the system

MODE

All modes

SYNTAX

INIT:<CRC16><CR>

REPLY

OKAY<CRC16><CR>

ERROR<ErrorCode><CRC16><CR>

PHSR

Return the number of assigned port handles and the port status for each one.

MODE

Setup

PREREQUISITE

INIT

SYNTAX

PHSR:<ReplyOption><CRC16><CR>

ReplyOption: 2 bytes. 01 for reporting port handles that need to be freed, and 02 for reporting port handles that are occupied but not initialized or enabled.

REPLY

<Numbers of PortHandle><nth PortHandle><nth

PortHandleStatus><CRC16><CR>

If no port handles, just REPLY 00<CRC16><CR>

Numbers of PortHandle: 2 bytes.

PortHandle: 2 bytes.

PortHandleStatus: 3 bytes.

PHRQ

Assign a port handle to a tool.

MODE

Setup

PREREQUISITE

INIT

SYNTAX

PHRQ:<HardwareDevice><SystemType><ToolType><PortNumber><Reserved><CRC16 ><CR>

HardwareDevice: 8 bytes.

SystemType: 1 byte.

ToolType: 1 byte. 0 for wired tool, and 1 for wireless, or * for wired.

PortNumber: 2 bytes. Reserved: 2 bytes.

REPLY

<PortHandle><CRC16><CR>

PortHandle: 2 bytes.

ERROR<ErrorCode><CRC16><CR>

PVWR

Assign a definition file to a tool.

MODE

Setup

PREREQUISITE

PHRQ or PHSR

SYNTAX

PVWR:<PortHandle><StartAddress><FileData><CRC16><CR>

PortHandle: 2 bytes. StartAddress: 4 bytes.

FileData: 128 bytes. When data is less than 128 bytes, pad '0' to

satisfy this condition.

A file larger than 128 bytes should be sent by calling this command for more than one time.

REPLY

OKAY<CRC16><CR>

ERROR<ErrorCode><CRC16><CR>

PINIT

Initialize a port handle

MODE

Setup

PREREQUISITE

PVWR or PHSR

SYNTAX

PINIT:<PortHandle><CRC16><CR>

PortHandle: 2 bytes.

REPLY

OKAY<CRC16><CR>

ERROR<ErrorCode><CRC16><CR>

PENA

Enables the reporting of transformation for a particular port handle.

MODE

Setup

PREREQUISITE

PINIT

SYNTAX

PENA:<PortHandle><TrackingPriority><CRC16><CR>

PortHandle: 2 bytes.

TrackingPriority: 1 byte. 'D' for Dynamic.

REPLY

OKAY<CRC16><CR>

ERROR<ErrorCode><CRC16><CR>

IRATE

Sets the illuminator rate.

MODE

Setup

PREREQUISITE

INIT

SYNTAX

IRATE:<IlluminatorRate><CRC16><CR>

IlluminatorRate: 1 byte. 0 for 20Hz, 1 for 30Hz, and 2 for 60Hz.

```
REPLY
OKAY<CRC16><CR>
ERROR<ErrorCode><CRC16><CR>
VSEL
Selects a characterized measurement volume.
MODE
Setup
PREREOUISITE
INIT
SYNTAX
VSEL:<VolumeNumber><CRC16><CR>
VolumeNumber: 1 byte.
REPLY
OKAY<CRC16><CR>
ERROR<ErrorCode><CRC16><CR>
SFLIST
Returns information about the supported features of the system.
MODE
Setup
SYNTAX
SFLIST:<ReplyOption><CRC16><CR>
ReplyOption: 2 byte.
   00 for summary of supported features,
   01 for number of active tool ports,
   02 for number of wireless tool ports,
   03 for characterized measurement volume,
   04 for number of wired tool ports, and
   05 for number of active wireless tool.
REPLY
<ReplyOptionData n><CRC16><CR>
ReplyOptionData 00: 8 bytes.
ReplyOptionData 01: 1 byte.
ReplyOptionData 02: 1 byte.
ReplyOptionData 03:
   <VolumeNumber><nth ShapeType><nth ShapeParameter><nth
WavelengthSupportedNumber><nth SupportedWavelength><LF>
   VolumeNumber: 1 byte.
   ShapeType: 1 byte.
   ShapeParameter: 10 parameters, 7 bytes each(a sign, 6 digits meaning
xxxx.xx).
   WavelengthSupportedNumber: 1 byte.
   SupportedWavelength: 1 byte per wavelength supported. 0 for 930nm, 1
```

```
for 880nm, and 4 for 870nm.
   <LF> is '\n'
ReplyOptionData 04: 1 byte.
ReplyOptionData 05: 1 byte.
ERROR<ErrorCode><CRC16><CR>
DSTART
Starts diagnostic mode
MODE
Setup
PREREQUISITE
INIT
SYNTAX
DSTART:<ReplyOption><CRC16><CR>
ReplyOption: optional, and 80 for resetting frame count to 0.
REPLY
OKAY<CRC16><CR>
ERROR<ErrorCode><CRC16><CR>
DSTOP
Stops diagnostic mode.
MODE
Diagnostic
PREREQUISITE
DSTART
SYNTAX
DSTOP: <CRC16><CR>
REPLY
OKAY<CRC16><CR>
ERROR<ErrorCode><CRC16><CR>
3D
Returns the latest 3D marker position of marker(s)
MODE
Diagnostic or Tracking
PREREQUISITE
IRED, only for active markers in Diagnostic mode.
SYNTAX
3D:<PortHandle><ReplyOption><CRC16><CR>
PortHandle: 2 bytes.
ReplyOption: 1 - 4 for single marker, 5 for up to 50 markers.
REPLY
<VisibleMarkerNumber><LF><ReplyOptionData n><CRC16><CR>
VisibleMarkerNumber: 3 bytes(a sign and 2 digits meaning xx).
```

```
ReplyOptionData 5: <Txn><Tyn><Tzn><LineSeparation><OutOfVolume><LF>
   Txn, Tyn, Tzn: 9 bytes each(a sign, 8 digits meaning xxxx.xxxx).
   LineSeparation: 4 bytes(a sign, 3 digits meaning x.xx).
   OutOfVolume: 1 byte. 0 for inside, and 1 for outside.
ERROR<ErrorCode><CRC16><CR>
TSTART
Starts Tracking mode.
MODE
Setup
PREREQUISITE
INIT
SYNTAX
TSTART:<ReplyOption><CRC16><CR>
ReplyOption: optional, and 80 for resetting frame count to 0.
REPLY
OKAY<CRC16><CR>
ERROR<ErrorCode><CRC16><CR>
TSTOP
Stops Tracking mode
MODE
Tracking
PREREQUISITE
TSTART
SYNTAX
TSTOP: <CRC16><CR>
REPLY
OKAY<CRC16><CR>
ERROR<ErrorCode><CRC16><CR>
TX
Returns the latest tool information.
MODE
Tracking
SYNTAX
TX:<ReplyOption><CRC16><CR>
ReplyOption: Optional. 4 bytes.
   0001 for transformation data by default,
   0002 for tool and marker information,
   0004 for 3D position of a single stray active marker,
   0008 for 3D positions of markers on tools,
   0800 for transformations not normally reported, and
   1000 for 3D positions of stray passive markers.
```

REPLY

<PortHandleNumber><PortHandle n><ReplyOptionData 0001>...<ReplyOptionData
0008><LF><ReplyOptionData 1000><SystemStatus><CRC16><CR>

PortHandleNumber: 2 bytes.

PortHandle: 2 bytes. ReplyOptionData 0001:

<Q0><Qx><Qy><Qz><Tx><Ty><Tz><Error><PortStatus><FrameNumber>

Q0, Qx, Qy, Qz: 6 bytes each(a sign, 5 digits meaning x.xxxx)

Tx, Ty, Tz: 7 bytes each(a sign, 6 digits meaning xxxx.xx)

Error: 6 bytes(a sign, 5 digits meaning x.xxxx).RMS error in mm.

PortStatus: 8 bytes. FrameNumber: 8 bytes.

Or MISSING<PortStatus><FrameNumber>

ReplyOptionData 0002 - 1000 could be null by default option.

SystemStatus: 4 bytes.

ERROR<ErrorCode><CRC16><CR>